

Datasheet

MYBPH monoclonal antibody (M02A), clone 1G18

Entrez GeneID: 4608**Gene Symbol:** MYBPH**Gene Alias:** -**Catalog Number:** H00004608-M02A**Regulation Status:** For research use only (RUO)**Product Description:** Mouse monoclonal antibody raised against a full-length recombinant MYBPH.**Clone Name:** 1G18**Immunogen:** MYBPH (AAH44226, 1 a.a. ~ 477 a.a)
full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.**Sequence:**

MMEKNTSEGPACSPREETASESAKVPTAEPPGEVAVSE
STREEQVPKPQAPAPQAPTASTATKPAPPSDEVPSAP
LLLTLDDVSSSSVTVSWEPPERLGRGLQGYVLELCR
EGASEWVPVSARPMVTQQTVRNALGDKFLLRVSA
VSSAGAGPPAMLDQPIHIRENIEAPKIRVPRHLRQTYIR
QVGETVNLQIPFQGKPKPQATWTHNGHALDSQRVSM
RTGDQDSILFIRSAQRSDSGRYELTVRVEDLEAKAVIDI
LVIEKPGPPSSIRLLDVWGCNAALQWTPPQDTGNTELL
GYMVQKADKKTGWFTVLERYPHTTCTISDLIIGNSYS
FRVFSENLGLSTSATVTKELAHQKADIAAKPKGFIER
DFSEAPSFTQPLADHTSTPGYSTQLFCSVRASPKPKII
WMKNKMEIQGNPKYRALSEQGVCTLEIRKPSPFDSGV
YTCKAINVLGEASVDCRLEVKASAAH

Host: Mouse**Reactivity:** Human**Applications:** ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Isotype:** IgM Kappa**Storage Buffer:** In ascites fluid**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.