

Datasheet

GPR139 polyclonal antibody

Gene Alias: GPRg1, PGR3

Catalog Number: PAB26294

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against synthetic peptide of GPR139.

Immunogen: A synthetic peptide corresponding to 17 amino acids at 1st cytoplasmic domain of human GPR139.

Host: Rabbit

Reactivity:

Bovine, Chicken, Dog, Gorilla, Hamster, Horse, Human, Monkey, Mouse, Pig, Rabbit, Rat

Applications: IHC-P

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: BLAST analysis of the peptide immunogen showed no homology with other human proteins, except GPR142 (65%).

Form: Liquid

Purification: Immunoaffinity chromatography

Recommend Usage: Immunohistochemistry

(Formalin/PFA-fixed paraffin-embedded sections) (6 ug/mL)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (0.1% sodium azide)

Storage Instruction: Store at 4°C. For long term storage store at -80°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 124274

Gene Symbol: GPR139