



Product Specification Sheet

VWR CRYOGENIC VIALS WITH EXTERNAL THREAD, LIP AND RED SILICONE O-RING

VWR International certifies that the specified products meet all VWR specifications and quality criteria for release. The products comply in every respect with the relevant particular specifications, drawings, relevant standards and regulations in force.

VWR International certifie que les produits mentionnés respectent toutes les spécifications et critères de qualité VWR. Les produits satisfont aux spécifications particulières, aux plans, aux normes et règlements en vigueur s'y rapportant.

VWR International bestätigt hiermit, daß die angegebenen Produkte alle VWR Spezifikationen und Qualitätskriterien zur Freigabe erfüllen. Die Produkte entsprechen in jeder Hinsicht den maßgeblichen Spezifikationen, den Zeichnungen sowie den relevanten Normen und gültigen Vorschriften.

Cat. No.	Capacity	Base	Ø ext.×H (mm)	Pk
479-1207	1.2 ml	Self standing	12.5×43	10 bags of 100 units
479-1209	2.0 ml	Round bottom	12.5×48	10 bags of 100 units
479-1208	2.0 ml	Self standing	12.5×49	10 bags of 100 units
479-1210	3.0 ml	Self standing	12.5×72	10 bags of 100 units
479-1217	4.0 ml	Self standing	12.5×76	10 bags of 100 units
479-1218	5.0 ml	Self standing	12.5×93	10 bags of 100 units

Material: **PP vial and closure, sterile**

Temperature range: **-196°C to +121°C Autoclavable at 121°C, for up to 30 minutes.**

Product compliance:

Resin conforms to FDA 21 CFR 177.1520

Resin conforms to USP Class VI

Plastics and colorants are in conformity with CONEG standards

Other Characteristics:

Certified RNase, DNase, Pyrogen and DNA free

Gamma radiation sterilized.

Printed graduations and writing area

Clear, no presence of contamination in plastic.

Leak proof testing in vacuum of 71.3cm/Hg.

Gas phase of Liquid nitrogen resistance testing.

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **August 2014**